



Bioprocessing

2D Labtainer Pro BioProcess Containers

A proven solution leading the industry in quality and reliability

2D Thermo Scientific™ Labtainer™ Pro BioProcess Containers (BPCs) are a proven solution that set the benchmark for industry excellence. With design flexibility throughout the portfolio, you can trust our high-quality films to help eliminate leaks and provide reliable liquid containment solutions from 50 mL to 20 L.

Quality:






- Proven, 11-quality film featuring leading single-use Thermo Scientific™ Aegis™ 5-14 film and Thermo Scientific™ CX5-14 film
- High level of assurance with 100% helium testing, automated manufacturing, and lot-based endotoxin and particulate testing
- Optimized drainage
- Minimized risk of leaks and simplified setup with the Thermo Scientific™ BioTitan™ Retention Device, which provides a 360° seal

Reliability:

- Reliable for product transportation and sampling, with scalability throughout the workflow
- High-quality products that are ready to ship
- Leading single-use manufacturing footprint for redundancy and supply assurance

Find it at eu.fishersci.com

Ordering information

Description	Cat. No	Fisher Scientific Cat. No
	Labtainer Pro BioProcess Container with BioTitan Retention, Aegis5-14 film, 3 ports, end-ported with AseptiQuik G connectors and body needleless Luer lock, 20 L	PL30024.14 17813956
	Labtainer Pro BioProcess Container with BioTitan Retention, Aegis5-14 film, 3 ports, end-ported with MPC body, MPC insert, and body needleless Luer lock, 20 L	PL30023.14 17853946
	Labtainer Pro BioProcess Container with BioTitan Retention, Aegis5-14 film, 2 ports, end-ported with insert and body Luer lock, 50 mL	PL30014.11 17823936
	Labtainer Pro BioProcess Container with BioTitan Retention, CX5-14 film, 3 ports, end-ported with MPC body, MPC insert, and body needleless Luer lock, 20 L	PL30010.14 17853926
	Labtainer Pro BioProcess Container with BioTitan Retention, Aegis5-14 film, 4 ports, end-ported with MPX, tri-clamps, and body needleless Luer lock, 10 L	PL30026.13 17893946

Learn more at eu.fishersci.com

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at **Belgium:** fishersci.be **Denmark:** fishersci.dk
Germany: fishersci.de **Ireland:** fishersci.ie **Italy:** fishersci.it
Finland: fishersci.fi **France:** fishersci.fr **Netherlands:** fishersci.nl
Norway: fishersci.no **Portugal:** fishersci.pt **Spain:** fishersci.es
Sweden: fishersci.se **Switzerland:** fishersci.ch **UK:** fishersci.co.uk

